

How to Install the CountMeds App for IOS



CountMeds
Compliance Strategies LLC for iOS.

1 thing you need to succeed

1. AppleID – you need your Apple ID and password.
 - Same that you use to download apps from the iTunes Store

4 Steps to help Us test CountMeds

1. Join iTunes Connect (Internal Testers Only)
2. Install Apple TestFlight
3. Install CountMeds
4. Set up CountMeds

1. Join iTunes Connect (Internal Testers Only)

If you are an Internal Tester you will receive an email like this from the iTunes Store inviting you to activate your iTunes Connect account.


If you are an External Tester, you can ignore this step

Internal vs External –how will you know?

- If you receive an email like this one then you are an Internal Tester.
- If you don't receive this email and only receive an email from TestFlight (see slide 6) then you are an External Tester


If you get this email then Tap “activate your account” in the email to start activating your iTunes Connect account

What you see on your phone




2. Install TestFlight

To install the TestFlight App you should go to the Apple App Store and search for "TestFlight"



Tap Install to complete this step




What you see on your phone



3. Install CountMeds

You will receive an email from TestFlight inviting you to test "CountMeds"

Tap the blue button labeled "View in TestFlight"



What you see on your phone



Install CountMeds

This is what the first screen of TestFlight will look like

Tap the blue button labeled "Continue"

What you see on your phone



18:14



Welcome to TestFlight



Accept Invites

Accept invites to beta programs and install the latest beta software.



Test Beta Apps

Test out the new functionality before it's released in the App Store.



Provide Feedback


Help developers build better apps by providing feedback and reporting bugs.

Continue

Apple ID: wm@coredevelopment.com. To use another Apple ID, go to Settings > iTunes & App Store.

Install CountMeds


Install the CountMeds App by selecting “Install”



What you see on your phone

Install CountMeds

Open CountMeds by selecting “OPEN”




What you see on your phone

4. Setup CountMeds

You must Setup the CountMeds application before you can help us test it.

Find the CountMeds Icon and tap on it.

What you see on your phone



Setup CountMeds

Step 1

You will come to a screen that allows you to start configuring the App.

Select "Start Here" to continue

**Your Count Meds app
is not configured.**

We need to set up your app before we can help you manage your prescription.

Start Here



Setup CountMeds

Step 2

Enter your Med-Check Patient ID.

It is 5 character code that can contain upper case letters, lower case letters and numbers such as

eD4Qh or
5PLw2

You should get your MedCheck ID from your clinic when you are first registered in the system

It is important to enter your MedCheck ID exactly!

Enter your Med-Check Patient ID to begin

Med-Check Patient ID

Next

Clinic

Please provide this ID to the Patient and record it into the Patient/Trial records.

Patient

Please provide your Med-Check Patient ID when you need assistance with this App. You should receive this from your medical provider.

What you see on your phone




Home

Setup CountMeds

Step 3

After you have entered your MedCheck ID,

Tap “Next” to continue to the next step



What you see on your phone

Setup CountMeds

Step 4

Next, enter your password and confirm your password by entering it a second time

Tap "Next" to continue to the next step

Please Log In

Med-Check Patient Id

pnTpT

Password

Login

What you see on your phone




Home

Setup CountMeds

Step 5

Select "Allow" when the CountMeds app asks you about sending notifications.

Notifications are necessary to send you reminders when it is time to photograph your pills



What you see on your phone


Setup CountMeds

Step 6

When you see this screen you are done configuring the application.

You should see your prescriptions that will be monitored with the CountMeds app

Congratulations!!!



What you see on your phone